

ARE MENT C.C.

Medicare-Hosquip

Medicare Hospital Equipment

Autoclaves,Ovens,Fridges and Washers.

ï¿½ï¿½ï¿½type=image,ypos=70,width=47.272222222,height=68.4388888889ï¿½ï¿½ï¿½q 134.00 0 0
194.00 28.35 619.54 cm /l2 Do Q

Fluid and blanket warmers and combination units made of stainless steel or mild steel epoxy coated. With digital electronic temperature control. With window or solid doors. Varying sizes

Autoclaves vertical or horizontal, varying sizes from 16lt cylindrical up to 400lt walk in. What is an autoclave - "An autoclave is a pressurized device designed to heat aqueous solutions above their boiling point at normal atmospheric pressure to achieve sterilization"

Vaccine refrigerators by Sibir Their cooling performance is unprecedented with efficient ice making capacity up to an ambient temperature of +50 C (+ 122 F). Sibir refrigerators have been tested according to the EPI standard E3/RF.2. The V170 is a powerful double door fridge offering high cooling and freezing performance. 185 litres gross capacity and a separate 47-litre freezer compartment.

[Fluid Warmers](#)

[Autoclaves](#)

[Vaccine Fridges](#)

[Surgeons Scrub Up
Sinks, Slop Hoppers](#)